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## INFORMATION MEDICALE IMPORTANTE EXRAY\*20, EXRAY\*30 et EXRAY\*35

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∅ DIAMÈTRE

## INSTRUCTION FOR USE

### IMPORTANT MEDICAL INFORMATION EXRAY\*20, EXRAY\*30 and EXRAY\*35

Reference : SUP\_7.016 version 2  
Date of last update : 11/2018

Date of first EC marking:  
-EXRAY\*20 : 2012  
-EXRAY\*30 / EXRAY\*35 : 2012

Product manufactured by :

**ADM**  
6 rue des Merisiers, 68920 Wettolsheim les Erlen – France  
Tel : +33 (0)3 89 23 75 64 // Fax : +33 (0)9 70 19 01 32  
[adv@synchromedical.com](mailto:adv@synchromedical.com)

**CAUTION : USA Federal law restricts this device to sale by or on the order of a physician.**

### INTENDED USE

The EXRAY\* screws system is intended for small bone extremity and large/long bone reconstruction fixation/or arthrodesis procedures for general use in skeletally mature individuals. The EXRAY\* screws system is intended to be permanently implanted without any other additional device and are delivered sterile.

### PURPOSE

EXRAY\* screws are designed to attach and stabilise bone segments in elective osteotomies on metatarsal bones and phalanges of the foot, in order to optimise the correct bone fusion of the segments concerned. These devices are only intended for use in the forefoot.

### DESCRIPTION:

The EXRAY\* screw system is intended for single use and is comprised of EXRAY\* 20, EXRAY\* 30 and EXRAY\* 35 screws.

EXRAY\* 20 screws comprise a single part implant in the form of a pin fitted with a threaded portion and a connector portion with a narrower diameter. EXRAY\*20 screws are non-compressive 2mm screws with a continuous thread.

EXRAY\* 30 screws comprise a single part implant. These 3mm self-compressive headless screws are available in several different lengths.

EXRAY\* 35 screws comprise a single part implant. These 3.5mm self-compressive headless screws are available in several different lengths.

EXRAY\* screws are manufactured in PEEK according to Standard ASTM F2026.

### INDICATIONS, CONTRAINDICATIONS AND POSSIBLE SIDE EFFECTS

**Indications :**  
The EXRAY\* screw system is indicated as an attachment system for small bone fractures or minor bone reconstruction in the following cases:

- monocortical or bicortical osteotomies of the foot or hand
- metatarsal or distal or proximal metacarpal osteotomies
- fusion of the first metatarsophalangeal and interphalangeal joint - attachment in the case of treatment of Hallux Valgus (e.g. scarf osteotomy, chevron osteotomy, etc.)
- Akin osteotomy

**CAUTION : to be used by or on the order of a licensed physician.**

TAKE NOTE : To be used on or according to the recommendation of a medical practitioner. This device may only be used by a doctor who has in-depth knowledge of the device, its intended use, and related instruments and surgical techniques.

### Contraindications:

- A non-exhaustive list of contraindications is as follows:
1. Any sign of generalised or local infection.
  2. Pathological obesity.
  3. Pregnancy.
  4. Any other medical or surgical condition that may compromise the success of surgery with instruments, such as the presence of malignant tumours, or serious congenital anomalies, an increase in sedimentation rates that cannot be attributed to other diseases, an increase in the number of white blood cells or a downward trend in such blood cells.
  5. Suspected or known allergy or intolerance to the implant's component materials.
  6. Any situation requiring the use of different materials
  7. Any case not listed in the indications.
  8. Any patient who is not willing to follow the postoperative instructions.
  9. Any patient in whom use of the implant may interfere with anatomical organs or some expected physiological function.
- The contraindications related to these devices are similar to those related to other osteosynthesis instruments. These osteosynthesis instruments have not been designed for, intended or sold for any use other than those indicated.

### Possible side effects:

In addition to risks associated with any surgery to the extremities, the list of possible undesirable side effects is as follows (non-exhaustive list):

1. Migration, fracture / rupture, expulsion, dislocation, or undesirable movement of one or more of the components
  2. Persistent pain and disability
  3. Bone loss related to short-circuiting of the constraints.
  4. The body's reaction to a foreign object owing to the presence of the implant (e.g. the development of a tumour, auto-immune disease and / or poor healing)
  5. Skin incompatibility with components where there is insufficient tissue coverage with pain and abnormal sensations due to the volume of the material.
  6. Superficial or deep infection with inflammatory reaction.
  7. Sepsis
  8. Pseudarthrosis or delay of bone fusion.
  9. Recurrence of the deformation, loss of correction
  10. Excessive loss of blood during the operation and / or haematoma.
  11. Lesions and / or vascular disorders (thrombosis)
  12. Peripheral circulatory complication (thrombosis and / or pulmonary embolism.
  13. Inability to resume the activities of normal everyday life.
  14. Compression of nearby tissues or organs.
- NOTE : Some additional surgical procedure may be required to remedy some of these possible undesirable events.

### INSERTION OF THE DEVICE

The EXRAY\* screw system must be used with instruments specially designed for the device according to the surgical technique concerned. It must not be used with any other instruments, unless this has specifically been recommended in some other ADM documentation, as using the system with other instruments poses the risk of incompatibility, and cannot be guaranteed. ADM shall not be held responsible in the case of using instruments that are not related to the implanted system.

**NEVER REUSE AN IMPLANT.**  
Reuse of an implant prevents the risk of contamination and loss of mechanical integrity.

### PACKAGING

The packaging of the device must be intact on receipt. Any device that has been damaged or any device whose packaging is damaged must not be used, and must rather be returned to ADM.

### STORAGE AND MAINTENANCE CONDITIONS

The implants must be stored in their original packaging, out of direct sunlight, humidity and extreme temperatures.

### MAGNETIC RESONANCE

The EXRAY\* screw system has not been evaluated for safety and compatibility in the MR environment. The EXRAY\* screw system has not been tested for heating or migration in the MR environment.

### CHECKING

The devices must always be checked before use. Those that present signs of damage or scratches on the surface must not be used. EXRAY\* implants are delivered sterile, and have been sterilised through gamma radiation. You need to carefully check that the round disc (sterilisation mark) on the packaging is red. Do not use an implant whose disc isn't red. You also need to check that the expiry date on the packaging has not passed. ADM disengages itself from any liability in the case of the use of implants past their expiry date. Any implant with damaged or torn packaging must not be re-sterilised but must rather be returned to ADM. TAKE NOTE : Do not re-sterilise an implant that has already been sterilised.

### WARNINGS AND PRECAUTIONS

Every surgery carries the risk of failure, as many external factors can compromise the outcome especially if the patient's state of health is already compromised. The implants must only be used in cases where bone fusion is deemed to be possible. No implant has been designed to withstand body weight in the long term without being accompanied by solid bone fusion. If a delay or lack of consolidation occurs, the implant may migrate and/or rupture. The following of preoperative and operative procedures, knowledge of surgical techniques, the selection of implants and their correct placement as well as appropriate selection and the patient's cooperation are factors leading to the operation's success. It has been shown that patients who smoke present a high risk of failure to obtain bone fusion. They need to be informed of this and any possible consequences. Obese patients, those who are malnourished and/or consume alcohol as well as patients with poor muscle or bone quality and/or those suffering from paralysis are also poor candidates for osteosynthesis surgery.

### OTHER PREOPERATIVE, PEROPERATIVE AND POSTOPERATIVE WARNINGS

**Implant selection**  
Selecting the suitable size of the implant for each patient is a crucial factor in the success of the operation, and it is important to refer to the surgical technique in this regard. Once implanted, the implant is subjected to repeated strain, and its strength is limited by adaptation of its geometrical structure to the size and shape of human bone. In order to minimise this strain, it is important to pay careful attention to the patient's selection criteria, the correct placement of the implant, and postoperative care. In fact, conversely, these constraints can lead to excessive strain placed on the material, resulting in deformation, rupture or loosening of the device likely to cause damage or the need to remove the implant prematurely. Any use in body areas not recommended by ADM cannot be guaranteed.

### Preoperative precautions

1. Only patients who meet the criteria described in the indications should be selected.
2. Patients who meet the criteria described in the contraindications should not be selected.
3. The implants must be handled and stored with the utmost care. They should not be scratched or damaged.
4. Before carrying out the procedure, the surgeon must check availability of the range of implants and associated instruments. They must personally handle all the devices before use in order to familiarise themselves with insertion of the implant.
5. Additional sterile components must be available in the event of unexpected need.

### Preoperative precautions

1. Instructions concerning the surgical technique should be read carefully. Any breakage or improper use of the implant may injure the patient or surgical team.
3. Specific instrumentation must be used to insert the device.
4. Take note: Do not use a device of the incorrect size (diameter and/or length), as this may damage the surrounding tissue.
6. Where possible and necessary, an imaging system may be used to check on placement of the device.

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